



Components of Lateral Flow



How Lateral Flow Tests Works

Test Line Capture goat antibody
against a specific agent

Control line Anti-species (anti-goat) antibody



Anthrax Test

False Positive Results

Test Line Capture goat antibody
against a specific agent

Control line Anti-species (anti-goat) antibody



Anthrax Test

Hook Effect or False Negative Results

Test Line Capture goat antibody
against a specific agent


Control line Anti-species (anti-goat) antibody



Anthrax Test

Understanding the Results

Positive



Two lines present in the result window
regardless of the intensity

Invalid



One line present at the Test line

Negative



One line present at the control line


Invalid



no lines present

Results may appear in as little as 3 minutes, with a maximum run time of 15

Understanding the Results



Positive

Two lines present in the result window.

Weak Positive

Two lines present in the result window regardless of the intensity

Negative

One line present at the control line

Invalid

One line present at the Test line

Invalid

no lines present

Results may appear in as little as 3 minutes, with a maximum run time of 15